

— ABLYBIO, Help Your Research



# Phospho-RIP3 (Ser164/Thr165) (YD34895) Rabbit mAb

货号: AYD16726

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	L929
细胞定位	Cytoplasm, cytosol, Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The product of this gene is a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases, and contains a C-terminal domain unique from other RIP family members. The encoded protein is predominantly localized to the cytoplasm, and can undergo nucleocytoplasmic shuttling dependent on novel nuclear localization and export signals. It is a component of the tumor necrosis factor (TNF) receptor-I signaling complex, and can induce apoptosis and weakly activate the NF-kappaB transcription factor.
基因ID	56532
基因名	Ripk3
Swiss	Q9QZL0 ( <a href="https://www.uniprot.org/uniprotkb/Q9QZL0/entry">https://www.uniprot.org/uniprotkb/Q9QZL0/entry</a> )
别名	Phospho-RIP3 (Ser164/Thr165) (YD34895),Phospho-RIP3 (Ser164/Thr165) (YD34895) Rabbit mAb,Ripk3,RIP-like protein kinase 3,Receptor-interacting protein 3,Rip3

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))