

CDKN2A/p16INK4a (YD33739) Rat mAb

货号: **AYD22628**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm, Nucleus
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Acts as a negative regulator of the proliferation of normal cells by interacting strongly with CDK4 and CD K6. This inhibits their ability to interact with cyclins D and to phosphorylate the retinoblastoma protein
基因ID	12578
基因名	Cdkn2a
Swiss	P51480
别名	CDKN2A/p16INK4a (YD33739)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn