

— ABLYBIO, Help Your Research



# Pentraxin 3 (YD34738) Rat mAb

货号: **AYD23627**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375,HL-60,HepG2,Jurkat,NIH/3T3,Mouse liver,Rat brain
细胞定位	Secreted
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the pentraxin protein family. The expression of this protein is induced by inflammatory cytokines in response to inflammatory stimuli in several mesenchymal and epithelial cell types, particularly endothelial cells and mononuclear phagocytes. The protein promotes fibrocyte differentiation and is involved in regulating inflammation and complement activation. It also plays a role in angiogenesis and tissue remodeling. The protein serves as a biomarker for several inflammatory conditions.
基因ID	5806
基因名	PTX3
Swiss	P26022 ( <a href="https://www.uniprot.org/uniprotkb/P26022/entry">https://www.uniprot.org/uniprotkb/P26022/entry</a> )
别名	Pentraxin 3 (YD34738),Pentraxin 3 (YD34738) Rat mAb,PTX3,Pentaxin-related protein PTX3,Tumor necrosis factor alpha-induced protein 5,Tumor necrosis factor-inducible gene 14 protein,TNFAIP5,TSG14

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))