

— ABLYBIO, Help Your Research



# beta II Tubulin (YD35288) Mouse mAb

货号: **AYD24177**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6,U-251MG,Mouse brain
细胞定位	Cytoplasm, cytoskeleton
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Microtubules, key participants in processes such as mitosis and intracellular transport, are composed of heterodimers of alpha- and beta-tubulins. The protein encoded by this gene is a beta-tubulin. Defects in this gene are associated with complex cortical dysplasia with other brain malformations-5. Two transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2015]
------	---

基因ID	347733
基因名	TUBB2B
Swiss	Q9BVA1 ( <a href="https://www.uniprot.org/uniprotkb/Q9BVA1/entry">https://www.uniprot.org/uniprotkb/Q9BVA1/entry</a> )
别名	beta II Tubulin (YD35288),beta II Tubulin (YD35288) Mouse mAb,TUBB2B

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)