

— ABLYBIO, Help Your Research



RHOB Rabbit mAb

货号: **AYM2006**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SK-BR-3,HeLa,mouse brain,mouse kidney,rat brain,rat lung
细胞定位	Cell membrane,Cleavage furrow,Late endosome membrane,Lipid-anchor,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human RHOB.
------	--

靶点信息

研究背景	RHOB (ras homolog family member B) is a protein-coding gene. GO annotations related to this gene include GDP binding and GTP binding. An important paralog of this gene is RND1.
基因ID	388
基因名	RHOB
Swiss	P62745
别名	RHOB,RHOB Rabbit mAb,Rho cDNA clone 6,ARH6,ARHB

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)