

— ABLYBIO, Help Your Research



# beta Tubulin Antibody

货号: **AYM2016**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>ICC:</b> 1:50 - 1:200 <b>FC:</b> 1:20 - 1:50
理论分子量	50kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,293T,SH-SY5Y,Mouse brain,Rat spleen
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human beta Tubulin.
------	--

## 靶点信息

研究背景	This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and 13. [provided by RefSeq, Jun 2014]
基因ID	203068
基因名	TUBB
Swiss	P07437
别名	CDCBM6,CSCSC1,M40,OK/SW-cl.56,TUBB1,TUBB5,beta Tubulin Antibody,TUBB,Tubulin beta-5 chain,beta Tubulin

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))