

pro Caspase 9 Rabbit mAb

货号: AYM28587

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	17kDa/30kDa/36kDa/46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,K-562,Mouse liver
细胞定位	apoptosome,cytoplasm,cytosol,mitochondrion,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human pro Caspase 9.
序列	VGALESLRGNADLAYILSMPCGHCLINNVNFCRESGLRTRTGSNIDCEKLRRRFSSLHFMVEVKGDLTAKKMWLALLELA QQDHGALDCCVVVILSHGCQASHLQFPGAVYGTGCPVSVEKIVNIFNGTSCPSLGGKPKLFFIQACGGEQKDHGFEVA STSPEDSPGSNPEPDATPFQEGLRFTDQLDAISSLPTPSDIFVSYSTFPGFVSWRDPKSGSWYVETLDDIFEQWAHSEDL QSLLLRVANAVSVKGIYKQMPGCFNFLRKKLFFKTS


靶点信息

研究背景	This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein can undergo autoproteolytic processing and activation by the apoptosome, a protein complex of cytochrome c and the apoptotic peptidase activating factor 1; this step is thought to be one of the earliest in the caspase activation cascade. This protein is thought to play a central role in apoptosis and to be a tumor suppressor. Alternative splicing results in multiple transcript variants.
基因ID	842
基因名	CASP9
Swiss	P55211
别名	CASP9;APAF-3;APAF3;ICE-LAP6;MCH6;PPP1R56;caspase-9;casp9;Caspase 9

产品验证



Immunohistochemical analysis of paraffin-embedded human colon, using pro Caspase 9 Antibody.



Western blot analysis of pro Caspase 9 expressed in Jurkat, K-562, Mouse liver using pro Caspase 9 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn