

— ABLYBIO, Help Your Research



LGR5/GPR49 Rabbit mAb

货号: **AYM28607**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP FC
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	100kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse ovary,Mouse brain,Rat testis,Rat brain
细胞定位	Cell membrane,Golgi apparatus,Multi-pass membrane protein,trans-Golgi network membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human LGR5/GPR49.
------	---

靶点信息

研究背景	The protein encoded by this gene is a leucine-rich repeat-containing receptor (LGR) and member of the G protein-coupled, 7-transmembrane receptor (GPCR) superfamily. The encoded protein is a receptor for R-spondins and is involved in the canonical Wnt signaling pathway. This protein plays a role in the formation and maintenance of adult intestinal stem cells during postembryonic development. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]
基因ID	8549
基因名	LGR5
Swiss	O75473 (https://www.uniprot.org/uniprotkb/O75473/entry)
别名	LGR5/GPR49,LGR5/GPR49 Rabbit mAb,LGR5,G-protein coupled receptor 49,G-protein coupled receptor 67,G-protein coupled receptor HG38,GPR49,GPR67

产品验证



Western blot analysis of LGR5/GPR49 expressed in Mouse ovary,Rat testis,Mouse brain using LGR5/GPR49 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)