

— ABLYBIO, Help Your Research



Integrin alpha V Rabbit mAb

货号: **AYM28609**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP ICC/IF
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	125kDa
实测分子量	125kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,MCF7,NIH/3T3,C6,Mouse brain,Mouse kidney,Rat kidney,Rat lung
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Integrin alpha V.
------	---


靶点信息

研究背景	The product of this gene belongs to the integrin alpha chain family. Integrins are heterodimeric integral membrane proteins composed of an alpha subunit and a beta subunit that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha V subunit. This subunit associates with beta 1, beta 3, beta 5, beta 6 and beta 8 subunits. The heterodimer consisting of alpha V and beta 3 subunits is also known as the vitronectin receptor. This integrin may regulate angiogenesis and cancer progression. Alternative splicing results in multiple transcript variants. Note that the integrin alpha 5 and integrin alpha V subunits are encoded by distinct genes. [provided by RefSeq, Oct 2015]
基因ID	3685
基因名	ITGAV
Swiss	P06756 (https://www.uniprot.org/uniprotkb/P06756/entry)
别名	Integrin alpha V,antigen identified by monoclonal L230,CD 51,CD51,DKFZp686A08142,integrin alpha V beta 3,integrin alpha-V,Integrin alpha-V light chain,integrin alphaVbeta3,integrin,alpha V (vitronectin receptor,alpha polypeptide,antigen CD51),ITAV_HUMAN,ITGAV,MSK 8,Msk8,Vitronectin receptor subunit alpha,V NRA,VTNR,Integrin alpha V Rabbit mAb,Vitronectin receptor

产品验证



Western blot analysis of Integrin alpha V expressed in Mouse brain,Rat kidney,A-549 using Integrin alpha V Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/protein s: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunofluorescence analysis of HeLa cells labelling Integrin alpha V with . The cells were fixed with cold 100% methanol (10min, 4°C) and blocked in 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween 20 for 1h. The cells were then incubate with Goat Anti-Rabbit IgG Fluor488 (1:200, shown in green) at room temperature for 1h.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)