

— ABLYBIO, Help Your Research



Syndecan 3 Rabbit mAb

货号: **AYM28671**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	45kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,MCF7,A-549,Mouse lung,Mouse spleen,Rat lung
细胞定位	cell surface,Golgi lumen,lysosomal lumen,microspike,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Syndecan 3.
------	---

靶点信息

研究背景	The protein encoded by this gene belongs to the syndecan proteoglycan family. It may play a role in the organization of cell shape by affecting the actin cytoskeleton, possibly by transferring signals from the cell surface in a sugar-dependent mechanism. Allelic variants of this gene have been associated with obesity.
------	---

基因ID	9672
基因名	SDC3
Swiss	O75056 (https://www.uniprot.org/uniprotkb/O75056/entry)
别名	Syndecan 3,Syndecan 3 Rabbit mAb,SDC3,KIAA0468

产品验证



Western blot analysis of Syndecan 3 expressed in Mouse lung,Rat lung,U-251MG using Syndecan 3 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)