

# Calmegin Rabbit mAb

货号: **AYM28681**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	<a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a> <a href="#">IP</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:20 - 1:50
理论分子量	100kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse testis,Rat testis
细胞定位	endoplasmic reticulum
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Calmegin.
序列	FGPDKCGEDYKLFHFRHKHPKTGVFEEKHAKPPDVDLKKFFTDKTHLYTLVMNPDDTFEVLVDQTVVNKGSLEDVVP PIKPPKEIEDPNDKKPEEWDERAKIPDPSAVKPEDWDESEPAQIEDSSVVKPAGWLDDEPKFIPDPNAEKPDWNETDGD EWEAPQILNPA

## 靶点信息

研究背景	Calmegin is a testis-specific endoplasmic reticulum chaperone protein. CLGN may play a role in spermatogenesis and infertility. [provided by RefSeq, Jul 2008]
基因ID	1047

基因名	CLGN
Swiss	O14967
别名	

## 产品验证



Western blot analysis of Calmegin expressed in PC-3, Mouse testis, Rat testis using Calmegin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)