

BMP2 Rabbit mAb

货号: **AYM28768**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse large intestine,Mouse stomach,Rat liver,HeLa,A-549,Rat spleen
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human BMP2.
序列	QAKHKQRKRLKSSCKRHPLYVDFSDVGVNDWIVAPPGYHAFYCHGECFPPLADHLNSTNHAIVQTLVNSVNSKIPKACCV PTELSAISMPLYLDENEKVVLLKNYQDMVVEGCGCR

靶点信息

研究背景	This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMA D family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which plays a role in bone and cartilage development. Duplication of a regulatory region downstream of this gene causes a form of brachydactyly characterized by a malformed index finger and second toe in human patients.
基因ID	650

基因名	BMP2
Swiss	P12643
别名	BDA2;BMP2A;BMP2

产品验证



Western blot analysis of BMP2 expressed in HT-29, Mouse large intestine, Mouse stomach using BMP2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn