

— ABLYBIO, Help Your Research



BIRC2 Rabbit mAb

货号: **AYM28769**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200
理论分子量	70kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,BxPC-3,HeLa,Jurkat
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human BIRC2 .
------	---

靶点信息

研究背景	The protein encoded by this gene is a member of a family of proteins that inhibits apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. This encoded protein inhibits apoptosis induced by serum deprivation and menadi one, a potent inducer of free radicals. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]
基因ID	329
基因名	BIRC2
Swiss	Q13490 (https://www.uniprot.org/uniprotkb/Q13490/entry)
别名	BIRC2,BIRC2 Rabbit mAb,Cellular inhibitor of apoptosis 1,IAP homolog B,Inhibitor of apoptosis protein 2,RING finger protein 48,RING-type E3 ubiquitin transferase BIRC2,TNFR2-TRAF-signaling complex protein 2,API1,MIHB,RNF48

产品验证



Western blot analysis of BIRC2 expressed in HepG2,293T,BxPC-3 using BIRC2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)