

Aggrecan Rabbit mAb

货号: **AYM28785**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	250kDa
实测分子量	kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Aggrecan.
序列	WPPTGAATEESTEGPSATEVPSASEEPSSEVPFPSEEPSSEEPFPSVRPFPSVELFPSEEPFPSKEPSPSEEPSASEEPTTP SPPVPSWTELPSSGEES

靶点信息

研究背景	This gene is a member of the aggrecan/versican proteoglycan family. The encoded protein is an integral part of the extracellular matrix in cartilagenous tissue and it withstands compression in cartilage. Mutations in this gene may be involved in skeletal dysplasia and spinal degeneration. Multiple alternatively spliced transcript variants that encode different protein isoforms have been observed in this gene. [provided by RefSeq, Jul 2008]
基因ID	176

基因名	ACAN
Swiss	P16112
别名	AGC1; AGCAN; CSPG1; CSPGCP; MSK16; SEDK

产品验证



Immunohistochemical analysis of paraffin-embedded human cartilage tissue, using Aggrecan A antibody.

实验步骤

访问官网浏览详情: www.ablybio.cn