

— ABLYBIO, Help Your Research



# Podoplanin Rabbit mAb

货号: **AYM28815**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	36kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,U-87MG,Mouse lung
细胞定位	Cell projection,Membrane,Single-pass type I membrane protein,filopodium membrane,lamellipodium membrane,microvillus membrane,ruffle membrane
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Podoplanin.
------	---

## 靶点信息

研究背景	This gene encodes a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor. The specific function of this protein has not been determined but it has been proposed as a marker of lung injury. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]
基因ID	10630
基因名	PDPN
Swiss	Q86YL7 ( <a href="https://www.uniprot.org/uniprotkb/Q86YL7/entry">https://www.uniprot.org/uniprotkb/Q86YL7/entry</a> )
别名	Podoplanin,Podoplanin Rabbit mAb,PDPN,Aggrus,Glycoprotein 36,PA2.26 antigen,T1-alpha,GP36

## 产品验证



Western blot analysis of Podoplanin expressed in U2OS,U-87MG,Mouse lung using Podoplanin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)