

Cytochrome C Rabbit mAb

货号: **AYM28941**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	11kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SH-SY5Y,BT-474
细胞定位	Mitochondrion intermembrane space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Cytochrome C.
序列	MGDVEKGKKIFIMKCSQCHTVEKGGKHKTGPNLHGLFGRKTGQAPGYSYAANKNKGIIWGEDTLMEYLENPKKYIPGTK MIFVGIKKKEERADLIAYLKKATNE

靶点信息


研究背景	This gene encodes a small heme protein that functions as a central component of the electron transport chain in mitochondria. The encoded protein associates with the inner membrane of the mitochondrion where it accepts electrons from cytochrome b and transfers them to the cytochrome oxidase complex. This protein is also involved in initiation of apoptosis. Mutations in this gene are associated with autosomal dominant nonsyndromic thrombocytopenia. Numerous processed pseudogenes of this gene are found throughout the human genome.
------	--

基因ID	54205
基因名	CYCS
Swiss	P99999
别名	CYCS;CYC;HCS;THC4;cytochrome c

产品验证



Immunohistochemical analysis of paraffin-embedded human colonic cancer, using Cytochrome C Antibody.



Western blot analysis of Cytochrome C expressed in 293T,SH-SY5Y,BT-474 using Cytochrome C Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates /proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn