

USP39 Rabbit mAb

货号: **AYM29007**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	65kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human USP39.
序列	KNPTIVNFPITNVDLREYLSEEVQAVHKNTTYDLIANIVHDGKPSSEGSYRIHVLHHGTGKWYELQDLQVTDILPQMILTSEAYI QIWKRRDNDQETNQGA

靶点信息

研究背景	Plays a role in pre-mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the pre-catalytic spliceosome. Regulates AURKB mRNA levels, and thereby plays a role in cytokinesis and in the spindle checkpoint. Does not have ubiquitin-specific peptidase activity.
基因ID	10713

基因名	USP39
Swiss	Q53GS9
别名	USP39; 65K; CGI-21; HSPC332; SAD1; SNRNP65; ubiquitin specific peptidase 39

产品验证



Western blot analysis of USP39 expressed in Mouse liver using USP39 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn