

CRABP1 Rabbit mAb

货号: **AYM29050**

产品信息

| | |
|-------|-----------------------------------------------------------------------------------------------|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IP |
| 推荐浓度 | WB: 1:500 - 1:2000 IP: 1:20 - 1:50 |
| 理论分子量 | 14kDa |
| 实测分子量 | 14kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | MCF7,NIH/3T3,Rat testis,Mouse testis |
| 细胞定位 | Cytoplasm |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|-----------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human CRABP1. |
|------|-----------------------------------------------------------|

靶点信息

| | |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a specific binding protein for a vitamin A family member and is thought to play an important role in retinoic acid-mediated differentiation and proliferation processes. It is structurally similar to the cellular retinol-binding proteins, but binds only retinoic acid at specific sites within the nucleus, which may contribute to vitamin A-directed differentiation in epithelial tissue. [provided by RefSeq, Jul 2008] |
| 基因ID | 1381 |
| 基因名 | CRABP1 |
| Swiss | P29762 |

| | |
|----|--------|
| 别名 | CRABP1 |
|----|--------|

产品验证



Western blot analysis of CRABP1 expressed in Mouse testis, MCF7, NIH/3T3 using CRABP1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn