

— ABLYBIO, Help Your Research



TWEAKR Rabbit mAb

货号: **AYM29121**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	14kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,PC-3
细胞定位	plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human TWEAKR.
------	---

靶点信息

研究背景	TNFRSF12A (TNF Receptor Superfamily Member 12A) is a Protein Coding gene. Diseases associated with TNFRSF12A include Multiple Sclerosis and Stroke, Ischemic. Among its related pathways are ERK Signaling and CDK-mediated phosphorylation and removal of Cdc6
基因ID	51330
基因名	TNFRSF12A
Swiss	Q9NP84 (https://www.uniprot.org/uniprotkb/Q9NP84/entry)
别名	TWEAKR,TWEAKR Rabbit mAb,TNFRSF12A,Fibroblast growth factor-inducible immediate-early response protein 14,Tweak-receptor, FN14

产品验证



Western blot analysis of TWEAKR expressed in HepG2,PC-3 using TWEAKR Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)