


基因ID	3849
基因名	KRT2
Swiss	P35908
别名	KRT2;CK-2e;K2e;KRT2A;KRT2E;KRTE;keratin 2

产品验证



Western blot analysis of Cytokeratin 2e expressed in HeLa, Mouse liver using Cytokeratin 2e Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn