

HDAC9 Rabbit mAb

货号: **AYM29234**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	57-65kDa/97-117kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RD,K-562
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human HDAC9.
序列	MHSMISSVDVKSEVPVGLPEISPLDLRTDLRMMMPVVDVPPVREKQLQQEQLLIQQQQIQKQLLIAEFQKQHENLTRQHQ AQLQEHIKELLAIKQQQELL


靶点信息

研究背景	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events . Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene has sequence homology to members of the histone deacetylase family. This gene is orthologous to the Xenopus and mouse MITR genes. The MITR protein lacks the histone deacetylase catalytic domain. It represses MEF2 activity through recruitment of multicomponent corepressor complexes that include CtBP and HDACs. This encoded protein may play a role in hematopoiesis. Multiple alternatively spliced transcripts have been described for this gene but the full-length nature of some of them has not been determined.
基因ID	9734
基因名	HDAC9
Swiss	Q9UKV0
别名	HDAC9;HD7;HD7b;HD9;HDAC;HDAC7;HDAC7B;HDAC9B;HDAC9FL;HDRP;MITR

产品验证



Immunohistochemical analysis of paraffin-embedded human brain, using HDAC9 Antibody.



Western blot analysis of HDAC9 expressed in RD,K-562 using HDAC9 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn