

MMP9 Rabbit mAb

货号: **AYM29253**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human MMP9.
序列	APRQRQSTLVLFPGDLRTNLTDRQLAEEYLYRYGYTRVAEMRGESKSLGPALLLLQQLSLPETGELDSATLKAMRTPRCG VPDLGRFQTFEGDLKWHHHNITYWIQNYSEDLPRVIDDAFARAFALWSAVTPLTFTTRVYSRDADIVIQFGVAEHGDGYP FDGKDGLLAHAFPPGPGIQQDAHFDDDELWSLKGKVVVPTFRGNADGAACHFPFIFEGRSYSACTTDGRSDGLPWCST TANYDTDDRFGFCPSERLYTQDGNADGKPCQFPFIFQGSYSACTTDGRSDGYRWCATTANYDRDKLFGFCPTRADST VMGGNSAGELCVFPFTFLGKEYSTCTSEGRGDGRLWCATTSNFDSDKKWGFCDQGYSLFLVAAHEFGHALGLDHSSV PEALMYPMYRFTEGPPLHKDDVNGIRHLYGPRPEPEPRPPTTTTTPQPTAPPTVCPTGPPTVHPSERPTAGPTGPPSAGPTGP PTAGPSTATTVPLSPVDDACNVNIFDAIAEIGNQLYLFKDGKYWRFSEGRGSRPQGPFLIADKWPALPRKLDVFEERLSK KLFFFSGRQVWVYTGASVLGPRRLDKLGLGADVAQVTGALRSGRGKMLLFSGRRLWRFDVKAQMVDPRSASEVDRMF PGVPLDTHDVFQYREKAYFCQDRFYWRVSSRSELNQVDQVGYVTYDILQCPED


靶点信息

研究背景	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. The enzyme encoded by this gene degrades type IV and V collagens. Studies in rhesus monkeys suggest that the enzyme is involved in IL-8-induced mobilization of hematopoietic progenitor cells from bone marrow, and murine studies suggest a role in tumor-associated tissue remodeling.
基因ID	4318
基因名	MMP9
Swiss	P14780
别名	MMP9;CLG4B;GELB;MANDP2;MMP-9

产品验证



Immunohistochemical analysis of paraffin-embedded human colon cancer, using MMP9 Antibody.



Western blot analysis of MMP9 expressed in Mouse brain using MMP9 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn