

— ABLYBIO, Help Your Research



# Phospho-Lyn (Y396) Rabbit mAb

货号: **AYM29284**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200 <b>FC:</b> 1:20 - 1:50
理论分子量	56kDa/58kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,K-562,Mouse lung,Rat lung
细胞定位	Cell membrane,Cytoplasm,Golgi apparatus,Lipid-anchor,Membrane,Nucleus,perinuclear region
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Phospho-Lyn (Y396).
------	---

## 靶点信息

研究背景	This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	4067
基因名	LYN
Swiss	P07948 ( <a href="https://www.uniprot.org/uniprotkb/P07948/entry">https://www.uniprot.org/uniprotkb/P07948/entry</a> )
别名	Phospho-Lyn (Y396),Phospho-Lyn (Y396) Rabbit mAb,LYN,Lck/Yes-related novel protein tyrosine kinase,V-y es-1 Yamaguchi sarcoma viral related oncogene homolog,p53Lyn,p56Lyn,JTK8

## 产品验证



Western blot analysis of Phospho-Lyn (Y396) expressed in Mouse lung,U-87MG,HeLa using Phospho-Lyn (Y396) Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))