

— ABLYBIO, Help Your Research



# Syntaxin Rabbit mAb

货号: **AYM29426**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IP:</b> 1:20 - 1:50
理论分子量	
实测分子量	33kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Y79 cells, fetal human brain tissue, rat brain tissue
细胞定位	Membrane ; Single-pass, Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Syntaxin.
------	---

## 靶点信息

研究背景	
------	--

基因ID	3134
基因名	STX1B
Swiss	P61266 ( <a href="https://www.uniprot.org/uniprotkb/P61266/entry">https://www.uniprot.org/uniprotkb/P61266/entry</a> )
别名	Syntaxin,Syntaxin Rabbit mAb,STX1B,Syntaxin-1B1,Syntaxin-1B2,STX1B1,STX1B2

## 产品验证



Western blot analysis of Syntaxin expressed in Y79 cells, fetal human brain tissue, rat brain tissue using Syntaxin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)