

— ABLYBIO, Help Your Research



Renilla Luciferase Rabbit mAb

货号: **AYM29444**

产品信息

反应	Renilla, Luciferase
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Recombinant protein
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Renilla Renilla Luciferase .
------	--

靶点信息

研究背景	
基因ID	

基因名	
Swiss	P27652 (https://www.uniprot.org/uniprotkb/P27652/entry)
别名	Renilla Luciferase, Renilla Luciferase Rabbit mAb, Renilla-luciferin 2-monooxygenase, Renilla-type luciferase

产品验证



Western blot analysis of Renilla Luciferase expressed in Recombinant protein using Renilla Luciferase Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)