

— ABLYBIO, Help Your Research



# PDK2 Rabbit mAb

货号: **AYM29562**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:20 - 1:50 <b>FC:</b> 1:20 - 1:50
理论分子量	46kDa
实测分子量	46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse heart,Rat heart
细胞定位	Mitochondrion matrix
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human PDK2.
------	---

## 靶点信息

研究背景	This gene encodes a member of the pyruvate dehydrogenase kinase family. The encoded protein phosphorylates pyruvate dehydrogenase, down-regulating the activity of the mitochondrial pyruvate dehydrogenase complex. Overexpression of this gene may play a role in both cancer and diabetes. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2010]
基因ID	5164
基因名	PDK2
Swiss	Q15119 ( <a href="https://www.uniprot.org/uniprotkb/Q15119/entry">https://www.uniprot.org/uniprotkb/Q15119/entry</a> )
别名	PDK2,PDK2 Rabbit mAb,Pyruvate dehydrogenase kinase isoform 2,PDHK2

## 产品验证



Western blot analysis of PDK2 expressed in Mouse brain, Mouse heart, Rat heart using PDK2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)