

TRK fused gene Rabbit mAb

货号: **AYM29723**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	30kDa/31kDa/43kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	cytoplasm,cytosol,endoplasmic reticulum exit site
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human TRK fused gene.
序列	MNGQLDLSGKLIKAQLGEDIRRIPIHNEDITYDELVLMMQRVFRGKLLSNDEVTIKYKDEDGDLITIFDSSDLSFAIQCSRIL KLTLFVNGQPRPLESSQVKYLRRELIELRNKVNRLDLSLEPPGEPGPSTNIPENDTVDGREEKSASDSSGKQSTQVMAASM SAFDPLKNQDEINKNVMSAFGLTDDQVSGPPSAPAEDRSGTPDSIASSSSAHPPGVQPQPPYTGAQTQAGQIE

靶点信息

研究背景	There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene.
------	--

基因ID	10342
基因名	TFG
Swiss	Q92734
别名	TFG;HMSNP;SPG57;TF6;TRKT3

产品验证



Western blot analysis of TRK fused gene expressed in HeLa using TRK fused gene Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn