

— ABLYBIO, Help Your Research



# COL11A1 Rabbit mAb

货号: **AYM29860**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	167kDa/176kDa/181kDa/182kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-204
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human COL11A1.
------	--

## 靶点信息

研究背景	This gene encodes one of the two alpha chains of type XI collagen, a minor fibrillar collagen. Type XI collagen is a heterotrimer but the third alpha chain is a post-translationally modified alpha 1 type II chain. Mutations in this gene are associated with type II Stickler syndrome and with Marshall syndrome. A single-nucleotide polymorphism in this gene is also associated with susceptibility to lumbar disc herniation. Multiple transcript variants have been identified for this gene.
基因ID	1301
基因名	COL11A1
Swiss	P12107 ( <a href="https://www.uniprot.org/uniprotkb/P12107/entry">https://www.uniprot.org/uniprotkb/P12107/entry</a> )
别名	COL11A1, COL11A1 Rabbit mAb, COLL6

## 产品验证



Western blot analysis of COL11A1 expressed in A-204 using COL11A1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)