

# PEX11B Rabbit mAb

货号: **AYM30010**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	28kDa
实测分子量	25kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse kidney,Mouse brain,Rat testis,Rat lung,Rat brain
细胞定位	cytosol,mitochondrion,nucleoplasm,peroxisome
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PEX11B.
序列	DNVLWAGKSG LAPRVDQEKWAQRSFRYYLFS LIMNLSRDAYEIRLLMEQESSACSRRLKSGGGVPGGSETGGLGGPG TPG

## 靶点信息

研究背景	The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene.
基因ID	8799

基因名	PEX11B
Swiss	O96011
别名	PEX14B;PEX11-BETA;PEX11B

## 产品验证



Western blot analysis of PEX11B expressed in HepG2, Mouse kidney, Mouse brain using PEX11B Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)