

— ABLYBIO, Help Your Research



Flotillin 1 Rabbit mAb

货号: **AYM30046**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC IP |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50 |
| 理论分子量 | 40kDa |
| 实测分子量 | 48kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,A375,Mouse liver |
| 细胞定位 | Secreted,extracellular matrix,extracellular space |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|---|
| 抗原信息 | Recombinant fusion protein corresponding to Human Flotillin 1 . |
|------|---|


靶点信息

| | |
|-------|--|
| 研究背景 | Stabilizes the aggregates of proteoglycan monomers with hyaluronic acid in the extracellular cartilage matrix. |
| 基因ID | 1404 |
| 基因名 | HAPLN1 |
| Swiss | P10915 |
| 别名 | Flotillin 1, Flotillin 1 Rabbit mAb, HAPLN1, Cartilage-linking protein 1, Proteoglycan link protein, CRTL1 |

产品验证



Western blot analysis of Flotillin 1 expressed in HeLa, A375, Mouse liver using Flotillin 1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunohistochemical analysis of paraffin-embedded human uterus, using Flotillin 1 Antibody .

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)