

PRP19 Rabbit mAb

货号: **AYM30071**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC ICC FC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 FC: 1:20 - 1:50 |
| 理论分子量 | 55kDa |
| 实测分子量 | 55kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | SW620 |
| 细胞定位 | Cytoplasm,Lipid droplet,Nucleus,cytoskeleton,nucleoplasm,spindle |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | Recombinant fusion protein corresponding to Human PRP19. |
| 序列 | MSLICISNEVPEHPCVSPVSNHVYERRLIEKYIAENGTDPINNQLSEEQOLIDIKVAHPIRPKPPSATSIPAILKALQDEWDAV MLHSFTLRQQLQTTTRQELSHALYQHDAACRVIA RLTKEVTAAREALATLKPQAGLIVPQAVPSSQPSVVGAGEPMDLGELV GMTPEIIQKLQDKATVLTTERKKRGKTVPEELVKPEELSKYRQVASHVGLHSAS |

靶点信息


| | |
|------|---|
| 研究背景 | PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).[supplied by OMIM, Sep 2008] |
| 基因ID | 27339 |

| | |
|-------|---|
| 基因名 | PRPF19 |
| Swiss | Q9UMS4 |
| 别名 | PRPF19;NMP200;PRP19;PSO4;SNEV;UBOX4;hPSO4 |

产品验证



Immunohistochemical analysis of paraffin-embedded human cervix cancer, using PRPF19 Antibody.



Western blot analysis of PRPF19 expressed in SW620 using PRPF19 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn