

— ABLYBIO, Help Your Research



# PRP19 Rabbit mAb

货号: **AYM30071**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>ICC:</b> 1:50 - 1:200 <b>FC:</b> 1:20 - 1:50
理论分子量	55kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620
细胞定位	Cytoplasm,Lipid droplet,Nucleus,cytoskeleton,nucleoplasm,spindle
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human PRP19.
------	--


## 靶点信息

研究背景	PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).[supplied by OMIM, Sep 2008]
基因ID	27339
基因名	PRPF19
Swiss	Q9UMS4
别名	PRP19,PRP19 Rabbit mAb,PRPF19,Nuclear matrix protein 200,PRP19/PSO4 homolog,RING-type E3 ubiquitin transferase PRP19,Senescence evasion factor,NMP200

## 产品验证



Western blot analysis of PRP19 expressed in SW620 using PRP19 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunohistochemical analysis of paraffin-embedded human cervix cancer, using PRP19 Antibody.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)