

— ABLYBIO, Help Your Research



# ENTPD5 Rabbit mAb

货号: **AYM30090**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver,Rat kidney
细胞定位	endoplasmic reticulum,extracellular region,extracellular space
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ENTPD5.
------	---

## 靶点信息

研究背景	The protein encoded by this gene is similar to E-type nucleotidases (NTPases)/ecto-ATPase/apyrases. NTP ases, such as CD39, mediate catabolism of extracellular nucleotides. ENTPD5 contains 4 apyrase-conserved regions which is characteristic of NTPases. [provided by RefSeq, Jan 2009]
------	---

基因ID	957
基因名	ENTPD5
Swiss	O75356
别名	ENTPD5,ENTPD5 Rabbit mAb,CD39 antigen-like 4,ER-UDPase,Ectonucleoside triphosphate diphosphohydr olase 5,Guanosine-diphosphatase ENTPD5,Inosine diphosphate phosphatase ENTPD5,Nucleoside diphosp hatase,Uridine-diphosphatase ENTPD5,CD39L4,PCPH

## 产品验证



Western blot analysis of ENTPD5 expressed in Mouse liver,Rat liver,Rat kidney using ENTPD5 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)