

# Uroplakin3a Rabbit mAb

货号: **AYM30101**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC FC
推荐浓度	<b>IHC:</b> 1:50 - 1:200 <b>FC:</b> 1:20 - 1:50
理论分子量	17kDa/30kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HL-60,Mouse brain,Mouse kidney,Rat kidney
细胞定位	Endoplasmic reticulum membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Uroplakin3a.
序列	VNLQPQLASVTFATNNPTLTTVALEKPLCMFDSKEALTGTHEVYLVLVDSAIRNASVQDSTNTPLGSTFLQTEGGRTGP YKAVAFDLIPCSDLPSLDAIGDVSKASQILNAYLVRVGANGTCLWDPNFQGLCNAPLSAATEYRFKYVLVNMSTGLVEDQT LWSDPIRTNQLTPYSTIDTWPGRSSGG

## 靶点信息

研究背景	This gene encodes a member of the uroplakin family, a group of transmembrane proteins that form complexes on the apical surface of the bladder epithelium. Mutations in this gene may be associated with renal adysplasia. Alternatively spliced transcript variants have been described.
基因ID	7380

基因名	UPK3A
Swiss	O75631
别名	UPK3A;UP3A;UPIII;UPIIIA;UPK3

## 产品验证



Immunohistochemical analysis of paraffin-embedded human bladder, using Uroplakin III Antibody.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)