

— ABLYBIO, Help Your Research



Uroplakin3a Rabbit mAb

货号: **AYM30101**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC FC
推荐浓度	IHC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	17kDa/30kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HL-60,Mouse brain,Mouse kidney,Rat kidney
细胞定位	Endoplasmic reticulum membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Uroplakin3a.
------	--

靶点信息

研究背景	This gene encodes a member of the uroplakin family, a group of transmembrane proteins that form complexes on the apical surface of the bladder epithelium. Mutations in this gene may be associated with renal dysplasia. Alternatively spliced transcript variants have been described.
基因ID	7380
基因名	UPK3A
Swiss	O75631 (https://www.uniprot.org/uniprotkb/O75631/entry)
别名	Uroplakin3a,Uroplakin3a Rabbit mAb,UPK3A,Uroplakin III,UPK3

产品验证



Immunohistochemical analysis of paraffin-embedded human bladder, using Uroplakin III Antibody.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)