



基因ID	358
基因名	AQP1
Swiss	P29972
别名	AQP1;AQP-CHIP;CHIP28;CO

产品验证



Immunohistochemical analysis of paraffin-embedded human breast, using AQP1 Antibody.



Western blot analysis of AQP1 expressed in Mouse liver, Mouse heart, Rat lung using AQP1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn