

— ABLYBIO, Help Your Research



STUB1 Rabbit mAb

货号: **AYM30227**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	27kDa/34kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29,OVCAR3,HT-1080,Mouse brain,Mouse kidney,Rat brain,Rat heart
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human STUB1.
------	--

靶点信息

研究背景	This gene encodes a protein containing tetratricopeptide repeat and a U-box that functions as a ubiquitin ligase/cochaperone. The encoded protein binds to and ubiquitinates shock cognate 71 kDa protein (Hspa 8) and DNA polymerase beta (Polb), among other targets. Mutations in this gene cause spinocerebellar ataxia, autosomal recessive 16. Alternative splicing results in multiple transcript variants. There is a pseudo gene for this gene on chromosome 2.
基因ID	10273
基因名	STUB1
Swiss	Q9UNE7
别名	STUB1,STUB1 Rabbit mAb,Antigen NY-CO-7,CLL-associated antigen KW-8,Carboxy terminus of Hsp70-interacting protein,RING-type E3 ubiquitin transferase CHIP,STIP1 homology and U box-containing protein 1,C HIP

产品验证



Western blot analysis of STUB1 expressed in Mouse brain,Rat brain,HeLa using STUB1 Rabbit mAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)