

— ABLYBIO, Help Your Research



Chromogranin B Rabbit mAb

货号: **AYM30291**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	125kDa
实测分子量	125kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,PC-12,BxPC-3,SH-SY5Y,Mouse brain
细胞定位	endoplasmic reticulum lumen,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Chromogranin B.
------	---

靶点信息

研究背景	This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides. [provided by RefSeq, Jan 2009]
------	--

基因ID	1114
基因名	CHGB
Swiss	P05060 (https://www.uniprot.org/uniprotkb/P05060/entry)
别名	Chromogranin B,Chromogranin B Rabbit mAb,CHGB,Chromogranin-B,Secretogranin I,SCG1

产品验证



Western blot analysis of Chromogranin B expressed in Mouse brain,293T,PC-12 using Chromogranin B Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30 ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)