

— ABLYBIO, Help Your Research



# FAM50A Rabbit mAb

货号: **AYM30398**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,PC-12,NIH/3T3,Mouse placenta
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human FAM50A.
------	---

## 靶点信息

研究背景	This gene belongs to the FAM50 family. The encoded protein is highly conserved in length and sequence across different species. It is a basic protein containing a nuclear localization signal, and may function as a DNA-binding protein or a transcriptional factor. [provided by RefSeq, Sep 2009]
基因ID	9130
基因名	FAM50A
Swiss	Q14320 ( <a href="https://www.uniprot.org/uniprotkb/Q14320/entry">https://www.uniprot.org/uniprotkb/Q14320/entry</a> )
别名	FAM50A,FAM50A Rabbit mAb,Protein HXC-26,Protein XAP-5,DXS9928E,HXC26,XAP5

## 产品验证



Western blot analysis of FAM50A expressed in Mouse placenta,U-251MG,PC-12 using FAM50A Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))