

mTOR Rabbit mAb

货号: **AYM30418**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	288kDa
实测分子量	289kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,C6
细胞定位	Cytoplasm,Cytoplasmic side,Endoplasmic reticulum membrane,Golgi apparatus membrane,Lysosome,Mitochondrion outer membrane,Nucleus,PML body,Peripheral membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human mTOR.
序列	MLGTGPAAATTAATTSSNVSVLQQFASGLKSRNEETRAKAAKELQHYYVTMELREMSQEESTRFYDQLNHHIFELVSSSDA NERKGGILAIASLIGVEGGNATRIGRFANYLRNLLPSNDPVMEMASKAIGRLAMAGDTFTAIEYVEFEVKRALEWLGADRN EGRRHAAVLVLRRELAISVPTFFFQQVQPFDFNIFVAVWDPKQAIREGAVAALRACLILTTQREPQEMQKQWYRHTFEEAEK GFDETLAKEKGMNRDDRIHGALLILNELVRISSMEGERLREEMEEITQQQLVHDKYC


靶点信息

研究背景	The protein encoded by this gene belongs to a family of phosphatidylinositol kinase-related kinases. These kinases mediate cellular responses to stresses such as DNA damage and nutrient deprivation. This protein acts as the target for the cell-cycle arrest and immunosuppressive effects of the FKBP12-rapamycin complex. The ANGPTL7 gene is located in an intron of this gene.
基因ID	2475
基因名	MTOR
Swiss	P42345
别名	FRAP;FRAP1;FRAP2;RAFT1;RAPT1;SKS;mTOR;MTOR

产品验证



Immunohistochemical analysis of paraffin-embedded rat liver, using mTOR Antibody.



Western blot analysis of mTOR expressed in K-562,C6 using mTOR Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn