

WNT2 Rabbit mAb

货号: **AYM30420**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 40kDa |
| 实测分子量 | 34kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Jurkat,NIH/3T3,Rat brain |
| 细胞定位 | Secreted,extracellular matrix,extracellular space |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein corresponding to Human WNT2. |
| 序列 | SWWYMRATGGSSRVMCDNVPGLVSSQRQLCHRHPDVMRAISQGVAEWTAECQHQRQHRWNCNTLDRDHSFLFGRV LLRSSRESAFVYAISSAGVVFAITRACSQGEVKSCSDPKKMGSADSKGIFDWGGCSDNIDYGIKFAFAFVDAKERKGG DARALMNLHNNRAGRKAVKRFLKQECKCHGVSGSCTLRTCWLAMADFRKTGDYLWRKYNGAIQVVMNQDGTGFTVAN ERFKKPTKNDLVYFENSPDYCIRDREAGSLGTAGRVCNLTSRGMDSCEVMCCGRGYDTSHVTRMTKCGCKFWCCAVR CQDCLEALDVHTCKAPKNADWTTAT |

靶点信息

| | |
|------|--|
| 研究背景 | This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene. |
|------|--|

| | |
|-------|-----------------|
| 基因ID | 7472 |
| 基因名 | WNT2 |
| Swiss | P09544 |
| 别名 | WNT2;INT1L1;IRP |

产品验证



Western blot analysis of WNT2 expressed in Jurkat,NIH/3T3,Rat brain using WNT2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn