

Plectin Rabbit mAb

货号: **AYM30476**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	531kDa
实测分子量	531kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	cytoplasm,cytosol,extracellular exosome,focal adhesion,intermediate filament cytoskeleton,mitochondrial outer membrane,perinuclear region of cytoplasm,plasma membrane,Z disc
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Plectin.
序列	

靶点信息

研究背景	
基因ID	5339
基因名	PLEC

Swiss	Q15149
别名	EBS1; EBSMD; EBSND; EBSO; EBSOG; EBSPA; HD1; LGMD2Q; PCN; PLEC1; PLEC1b; PLTN

产品验证



Western blot analysis of Plectin expressed in HUVEC, mouse heart, rat lung using Plectin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn