

— ABLYBIO, Help Your Research



Plectin Rabbit mAb

货号: **AYM30476**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	531kDa
实测分子量	531kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	cytoplasm,cytosol,extracellular exosome,focal adhesion,intermediate filament cytoskeleton,mitochondrial outer membrane,perinuclear region of cytoplasm,plasma membrane,Z disc
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Plectin.
------	--

靶点信息

研究背景	
------	--

基因ID	5339
基因名	PLEC
Swiss	Q15149
别名	Plectin,Plectin Rabbit mAb,PLEC,Hemidesmosomal protein 1,Plectin-1,PLEC1

产品验证



Western blot analysis of Plectin expressed in HUVEC,mouse heart, rat lung using Plectin Rabbit mAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)