

Histone H2A.Z Rabbit mAb

货号: **AYM30485**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	13kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,MCF7,22Rv1,PC-12,HepG2,Mouse testis,Mouse thymus,293T
细胞定位	Chromosome,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Histone H2A.Z.
序列	MAGGKAGKDSGKAKTKAVSRSQRAGLQFPVGRHRHLKSRRTTSHGRVGATAAVYSAAILEYLTAEVLELAGNASKDLKVK RITPRHLQLAIRGDEELDSLKATIAGGGVIPHIHKSIGKKGQKTV


靶点信息

研究背景	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent member of the histone H2A family that is distinct from other members of the family. Studies in mice have shown that this particular histone is required for embryonic development and indicate that lack of functional histone H2A leads to embryonic lethality.
基因ID	3015
基因名	H2AFZ
Swiss	POC0S5
别名	H2AFZ;H2A.Z-1;H2A.z;H2A/z;H2AZ;histone H2A.Z

产品验证



Immunohistochemical analysis of paraffin-embedded human kidney, using Histone H2A.Z Antibody.



Western blot analysis of Histone H2A.Z expressed in Mouse testis, HL-60, MCF7 using Histone H2A.Z Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn