

Calmegin Rabbit mAb

货号: **AYM30496**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	100kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse testis,Rat testis
细胞定位	endoplasmic reticulum
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Calmegin.
------	---

靶点信息

研究背景	Calmegin is a testis-specific endoplasmic reticulum chaperone protein. CLGN may play a role in spermatogenesis and infertility. [provided by RefSeq, Jul 2008]
基因ID	1047
基因名	CLGN

Swiss	O14967
别名	Calmegin

产品验证



Western blot analysis of Calmegin expressed in PC-3, Mouse testis, Rat testis using Calmegin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn