

— ABLYBIO, Help Your Research



UCHL5IP Rabbit mAb

货号: **AYM30584**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Mouse testis,Rat testis
细胞定位	centrosome,cytosol,nucleolus,plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human UCHL5IP.
------	--

靶点信息

研究背景	This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012]
基因ID	55559
基因名	HAUS7
Swiss	Q99871 (https://www.uniprot.org/uniprotkb/Q99871/entry)
别名	UHL5IP,UHL5IP Rabbit mAb,HAUS7,26S proteasome-associated UCH37-interacting protein 1,UHL5-interacting protein,X-linked protein STS1769,UIP1

产品验证



Western blot analysis of UHL5IP expressed in HeLa, Jurkat, Mouse testis using UHL5IP Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)