

GATA2/3 Rabbit mAb

货号: **AYM30626**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	51kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse lung,Rat lung
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GATA2/3.
序列	MEVAPEQPRWMAHPAVLNAQHPDSSHHPGLAHNYMEPAQLLPDEVDVFFNHLD SQGNPYANPAHARARVSYSPA HAR LTGGQMCRPHLLHSPGLPWLDGGKAALSAAAAHHHPWTVSPFSKTPHPSAAGGPGGPLSVYPGAGGGSGGGSGSS VASLTPTAAHSGSHLFGFPPTPPKEVSPDPSTTGAASSPASSAGGSAARGEDKDG

靶点信息


研究背景	This gene encodes a member of the GATA family of zinc-finger transcription factors that are named for the consensus nucleotide sequence they bind in the promoter regions of target genes. The encoded protein plays an essential role in regulating transcription of genes involved in the development and proliferation of hematopoietic and endocrine cell lineages. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Mar 2009]
------	---

基因ID	2624,2625
基因名	GATA2,GATA3
Swiss	P23769/P23771
别名	DCML; IMD21; MONOMAC; NFE1B

产品验证



Immunohistochemical analysis of paraffin-embedded mouse kidney, using GATA2/3 Antibody.



Western blot analysis of GATA2/3 expressed in HeLa, Mouse lung, Rat lung using GATA2/3 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn