

# Cyclin H/p34 Rabbit mAb

货号: **AYM30734**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	37kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa cells, LNCaP cells, HT-1080 cells, NCCIT cells, MCF-7 cells, HSC-T6 cells, Jurkat cells, HepG2 cells, H EK-293 cells, 4T1 cells, NIH/3T3 cells, Caco-2 cells, A431 cells
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Cyclin H/p34.
序列	MYHNSSQKRHWTFSSSEQLARLRADANRKFRCCKAVANGKVLNPNPVPFLEPHEEMTLCKYYEKRLLEFCVFKPAMPRSVV GTACMYFKRFYLNNSVMEYH

## 靶点信息

研究背景	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with CDK7 kinase and ring finger protein MAT1. The kinase complex is able to phosphorylate CDK2 and CDC2 kinases, thus functions as a CDK-activating kinase (CAK). This cyclin and its kinase partner are components of TFIIF, as well as RNA polymerase II protein complexes. They participate in two different transcriptional regulation processes, suggesting an important link between basal transcription control and the cell cycle machinery. A pseudogene of this gene is found on chromosome 4. Alternate splicing results in multiple transcript variants .
基因ID	902
基因名	CCNH
Swiss	P51946
别名	CCNH; CAK; Cych; p34; p37; cyclin-H

## 产品验证



Western blot analysis of Cyclin H/p34 expressed in HeLa cells, LNCaP cells, HT-1080 cells using Cyclin H/p34 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates /proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)