

# PRPF4 Rabbit mAb

货号: **AYM30767**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	58kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse spleen,Mouse thymus,Rat liver
细胞定位	Nucleus speckle
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PRPF4.
序列	

## 靶点信息

研究背景	The protein encoded by this gene is part of a heteromeric complex that binds U4, U5, and U6 small nuclear RNAs and is involved in pre-mRNA splicing. The encoded protein also is a mitotic checkpoint protein and a regulator of chemoresistance in human ovarian cancer. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	9128

基因名	PRPF4
Swiss	O43172
别名	PRPF4;HPRP4;HPRP4P;PRP4;Prp4p;RP70;SNRNP60

## 产品验证



Western blot analysis of PRPF4 expressed in Mouse spleen,Rat liver,HeLa using PRPF4 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)