

— ABLYBIO, Help Your Research



PRPF4 Rabbit mAb

货号: **AYM30767**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	58kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse spleen,Mouse thymus,Rat liver
细胞定位	Nucleus speckle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human PRPF4.
------	--

靶点信息

研究背景	The protein encoded by this gene is part of a heteromeric complex that binds U4, U5, and U6 small nuclear RNAs and is involved in pre-mRNA splicing. The encoded protein also is a mitotic checkpoint protein and a regulator of chemoresistance in human ovarian cancer. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	9128
基因名	PRPF4
Swiss	O43172 (https://www.uniprot.org/uniprotkb/O43172/entry)
别名	PRPF4,PRPF4 Rabbit mAb,PRP4 homolog,U4/U6 snRNP 60 kDa protein,WD splicing factor Prp4,PRP4

产品验证



Western blot analysis of PRPF4 expressed in Mouse spleen,Rat liver,HeLa using PRPF4 Rabbit mAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)